



**UNITED STATES DEPARTMENT OF COMMERCE**  
**National Institute of Standards and Technology**  
Gaithersburg, Maryland 20899-0001

DATE: 01 December 2015

**Product Identifier**

**SRM Number:** 1508a

**SRM Name:** Benzoylcegonine (Cocaine Metabolite) in Freeze-Dried Urine

Under the U.S. Department of Labor, Occupational Safety and Health Administration (OSHA) 29 CFR 1910.1200, this Standard Reference Material (SRM) is NOT classified as a physical hazard or a health hazard, a simple asphyxiant, combustible dust, pyrophoric gas, or hazard not otherwise classified. There are no hazard pictograms, hazard statements or signal word associated with it. Safety Data Sheet information is not required. This document may be used in conjunction with your hazard communication program.

**Exemption:** 1910.1200(b)(6)(xii). This SRM is a biological material and should be considered a potential biological hazard.

**Description:** This SRM 1508a is intended primarily for validating methods used for the determination of benzoylcegonine (cocaine metabolite) in human urine. SRM 1508a consists of four bottles of freeze-dried urine: one bottle each of three different levels of benzoylcegonine plus one bottle of blank freeze-dried urine.

**Additional Notes for Biomaterials:** SRM 1508a IS INTENDED FOR RESEARCH USE. THIS IS A HUMAN-SOURCE MATERIAL. HANDLE PRODUCT AS A BIOHAZARDOUS MATERIAL CAPABLE OF TRANSMITTING INFECTIOUS DISEASE. Accordingly, this human urine based product should be handled at the Biosafety Level 2 or higher as recommended for any POTENTIALLY INFECTIOUS HUMAN SPECIMEN in the Centers for Disease Control and Prevention/National Institutes of Health Manual.

**Disposal:** SRM 1508a components and derived solutions should be disposed of in accordance with local, state, and federal regulations.

**Transport Information:** This material is not regulated by the U.S. Department of Transportation (DOT) and/or International Air Transportation Association (IATA).

**Disclaimer:** This document was prepared carefully, using current references. Users of this SRM should ensure that this document and the corresponding Certificate of Analysis in their possession are current. This can be accomplished by contacting the SRM Program: telephone (301) 975-2200; fax (301) 948-3730; e-mail [srmmsds@nist.gov](mailto:srmmsds@nist.gov); or via the Internet at <http://www.nist.gov/srm>.