

301.3 - Electrophoretic Mobility, E (suspension form)

SRM 1980 is intended for use in the calibration and evaluation of equipment used to measure electrophoretic mobility. It consists of a goethite suspension saturated with phosphate in a sodium perchlorate electrolyte solution.

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

SRM	1980	1992	1993
Description	Positive Electrophoretic Mobility (+ μ_E) Standard	Zeta Potential - Colloidal Silica (Nominal Mass Fraction 0.15 %)	Zeta Potential - Colloidal Silica (Nominal Mass Fraction 2.2 %)
Unit Size	(40 mL)	(4 x 5 mL)	(2 x 25 mL)
Property Certified	+ μ_E , 2.53 $\mu\text{m} \cdot \text{cm/V} \cdot \text{s}$	Mean electrophoretic mobility, -4.5 x 10 ⁻⁸ m ² V ⁻¹ s ⁻¹ Mean zeta potential, -58 mV	Mean electrophoretic mobility, -4.3 x 10 ⁻⁸ m ² V ⁻¹ s ⁻¹ Mean zeta potential, -56 mV

- Certified values are normal font
- Reference values are italicized
- Values in parentheses are for information only