

105.1 - Clinical Laboratory Materials (gas, liquid, and solid forms)

The following SRMs are for calibrating apparatus and validating analytical methods used in clinical and pathology laboratories. Additional information on the serum materials is given in table 105.2.

For further information see: [SP 260-36](#), [SP 260-72](#) and [SP 260-83](#)

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

SRM Description	900a	909c	911c	912a	913b	914a	915b	916a	917c	918c	919b	920	921	924a	927e	928	929a	937	955c	956d
Unit Size	Antiepilepsy Drugs in Frozen Human Serum (4 vials)	Frozen Human Serum (3 x 2 mL)	Cholesterol (2 g)	Urea (25 g)	Uric Acid (10 g)	Creatinine (10 g)	Calcium Carbonate (20 g)	Bilirubin (100 mg)	D-Glucose (Dextrose) (50 g)	Potassium Chloride General and Ion Activity Standard (Dried at 110 °C) (30 g)	Sodium Chloride (30 g)	D-Mannitol (50 g)	Cortisol (Hydrocortisone) (1 g)	Lithium Carbonate (30 g)	Bovine Serum Albumin (7% solution) (Total Protein Standard) (10 x 2.2mL)	Lead Nitrate (30 g)	Magnesium Gluconate (5 g)	Iron Metal (50 g)	Toxic Metals in Caprine Blood (4 vials)	Electrolytes in Frozen Human Serum (6 x 2.0 mL)
Purity/Constituent (mass fraction in %)	2 levels: phenobarbital, phenytoin, lamotrigine, topiramate	Cholesterol, Elements, Creatinine, Total Glycerides, Urea, Uric Acid, Protein	99.2	99.9	99.8	99.7	WCaCO ₃ 99.907 WCa 40.0104 WCO ₃ 59.923	98.3	99.7	WKCl 99.945 WK 52.421 WCl 47.5317	WNaCl 99.835 WCl ⁻ 60.564 WNa ⁺ 39.2747	99.8	98.9	99.867*	BSA Conc. 67.38 g/L	100.00	Mg Conc. 5.362	99.90	4 Levels: Pb, As, Cd, Hg 1 Level: ethyl- methyl- and inorganic Hg species	3 Levels: Ca, Cl, Li, Mg, K, Na 3 Levels: ionized calcium, Phosphate, and Phosphorus

* Conforms to National Committee for Clinical Laboratory Standards (NCCLS) specification ACC-1.

- Certified values are normal font
- Reference values are italicized
- Values in parentheses are for information only

105.1 - Clinical Laboratory Materials (gas, liquid, and solid forms)

The following SRMs are for calibrating apparatus and validating analytical methods used in clinical and pathology laboratories. Additional information on the serum materials is given in table 105.2.

For further information see: [SP 260-36](#), [SP 260-72](#) and [SP 260-83](#)

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

* Conforms to National Committee for Clinical Laboratory Standards (NCCLS) specification ACC-1.

- Certified values are normal font
- Reference values are italicized
- Values in parentheses are for information only

105.1 - Clinical Laboratory Materials (gas, liquid, and solid forms)

The following SRMs are for calibrating apparatus and validating analytical methods used in clinical and pathology laboratories. Additional information on the serum materials is given in table 105.2.

For further information see: [SP 260-36](#), [SP 260-72](#) and [SP 260-83](#)

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

SRM	965b	967a	968e	971	972a	998	1400	1486	1595	1598a	1950	1951c	1952a	2366a
Description	Glucose in Frozen Human Serum	Creatinine in Frozen Human Serum	Fat-Soluble Vitamins, Carotenoids, and Cholesterol in Human Serum	Hormones in Frozen Human Serum	Vitamin D Metabolites in Frozen Human Serum	Angiotensin I (Human)	Bone Ash	Bone Meal	Tripalmitin	Inorganic Constituents in Animal Serum	Metabolites in Frozen Human Plasma	Lipids in Frozen Human Serum	Cholesterol in Freeze-Dried Human Serum	Cytomegalovirus DNA (Townsend 147 BAC) for DNA Measurements
Unit Size	(set (8) (2 each conc))	(set(4) (2 each conc))	(set (3) (1 each conc.))	(2 x 5 mL)	(4 x 1 mL)	(0.5 mg)	(50 g)	(50 g)	(2 g)	(2 vials)	(5 x 1 mL)	(4 vials (2 each conc))	(set (6) (2 each conc))	(1 vial)
Purity/Constituent (mass fraction in %)	4 Levels: glucose	2 Levels: creatinine	3 Levels: vitamins, carotenoids, and cholesterol	cortisol, testosterone, progesterone, total thyroxine (T4), total 3,3',5-triiodothyronine (T3)	4 Levels: 25-hydroxyvitamin D ₂ , 25-hydroxyvitamin D ₃ , 3-epi-25-hydroxyvitamin D ₃	94.1	8 elements	8 elements	99.5	12 elements 5 elements	Cholesterol, Creatinine, Urea,Uric Acid, Homocysteine, Glucose, Total Glycides, 3 Hormones, 8 Fatty Acids, 12 Amino Acids, 9 Vitamins, 4 Elements, 2 carotenoids, Bilirubin, Selenium Species, Total protein, 6 PFCs, 19 Fatty Acids, 4 Amino Acids, 4 Vitamins, 3 Elements, 2 carotenoids,	2 Levels: Total Cholesterol, Total Glycerides	3 Levels: cholesterol	1 Level: cytomegalovirus (copies per micro liter)

* Conforms to National Committee for Clinical Laboratory Standards (NCCLS) specification ACC-1.

- Certified values are normal font
- Reference values are italicized
- Values in parentheses are for information only

105.1 - Clinical Laboratory Materials (gas, liquid, and solid forms)

The following SRMs are for calibrating apparatus and validating analytical methods used in clinical and pathology laboratories. Additional information on the serum materials is given in table 105.2.

For further information see: [SP 260-36](#), [SP 260-72](#) and [SP 260-83](#)

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

2378	2389a	2668	2669	2670a	2921
Fatty Acids in Frozen Human Serum (3 vials, 1 level each)	Amino Acids in 0.1 mol/L Hydrochloric Acid (5 x 1.2 mL)	Toxic Elements in Frozen Human Urine (Set (10) (5 each conc))	Arsenic Species in Frozen Human Urine (Set (10) (5 each conc.))	Toxic Elements in Freeze-Dried Urine (set (4) (2 each conc))	Human Cardiac Troponin Complex (5 x 115 uL)
3 Levels, Fatty acids	17 Amino Acids	2 Levels: 14 Elements 2 Levels: 9 Elements	2 levels: Arsenic species 2 levels: Total Arsenic	2 Levels: 14 Elements 10 Elements	Cardiac Troponin, cTnI, cTnT, cTnC

* Conforms to National Committee for Clinical Laboratory Standards (NCCLS) specification ACC-1.

- Certified values are normal font
- Reference values are italicized
- Values in parentheses are for information only

105.1 - Clinical Laboratory Materials (gas, liquid, and solid forms)

The following SRMs are for calibrating apparatus and validating analytical methods used in clinical and pathology laboratories. Additional information on the serum materials is given in table 105.2.

For further information see: [SP 260-36](#), [SP 260-72](#) and [SP 260-83](#)

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

SRM	2972a	2973	3667	3668	3669	3671	3672	3673	3950	8321	8323	8327	8642a	8671
Description	25-Hydroxyvitamin D Calibration Solutions	Vitamin D Metabolites in Frozen Human Serum (High Level)	Creatinine in Frozen Human Urine	Mercury, Perchlorate, and Iodide in Frozen Human Urine	Arsenic Species in Frozen Human Urine (Elevated Levels)	Nicotine Metabolites in Human Urine (Frozen)	Organic Contaminants in Smokers' Urine (Frozen)	Organic Contaminants in Non-Smokers' Urine (Frozen)	Vitamin B ₆ in Frozen Human Serum	Peptide Mixture of Proteomics	Yeast Protein Extract	Peptide Reference Material for Molecular Mass and Purity Measurements	FDA Saxitoxin Dihydrochloride Solution	NISTmAb, Humanized IgG1k Monoclonal Antibody
Unit Size	(20 x 1.2 mL)	(one vial)	(1 x 10 mL)	(Set (10) (5 each conc))	(5 x 1.5mL)	(3 x 10 mL)	(5 x 10 mL)	(5 x 10 mL)	(2 vials, 1 level each)	(3 vials)	(3 x 0.2 mL)	(set 3 (1 mg each))	(5 ampoules x 1.2 mL)	(800 microliters)
Purity/Constituent (mass fraction in %)	1 Level: 25-hydroxyvitamin D ₂ 2 Levels: 25-hydroxyvitamin D ₃ 1 Level: 3-Epi-25-hydroxyvitamin D ₃	25-hydroxyvitamin D ₃ 25-hydroxyvitamin D ₂ 3-epi-25-hydroxyvitamin D ₃ 24R ₂ 25-dihydroxyvitamin D ₂	61.8 mg/dL	2 Levels: Hg, Nitrate, Perchlorate 2 Levels: I, Thiocyanate	Arsenic Species	2 Levels: Nicotine, Cotinine 1 Level: 3-hydroxycotinine Other Nicotine Metabolites	Hydroxylated PAHs Phthalate Metabolites Phenol Metabolites VOC Metabolites	Hydroxylated PAHs Phthalate Metabolites Phenol Metabolites VOC Metabolites	2 Levels: Vitamin B ₆	440 Synthetic Peptides	0.191 mg/mL	Molecular Mass, Purity of Peptides, Amino Acid Sequence, for 3 synthetic peptides	0.0103	Mass Concentration (1 level) and Physicochemical Values

* Conforms to National Committee for Clinical Laboratory Standards (NCCLS) specification ACC-1.

- Certified values are normal font
- Reference values are italicized
- Values in parentheses are for information only