

105.2 - Serum and Plasma Materials (frozen, liquid, and lyophilized forms)

These SRMs serve a variety of clinical measurement needs. SRM [909c](#) is a frozen human serum for use in determining specified constituents. SRM [927d](#) is a bovine serum albumin in a sterile 7% solution for use in the calibration and standardization of procedures to analyze total serum protein. SRM [956c](#) is a frozen human serum for use in the calibration and standardization of procedures for the determination of specific electrolytes in either diluted or undiluted human serum or plasma. SRM [965b](#) is a frozen human serum for evaluating the accuracy of procedures used to determine glucose in human serum and to validate secondary reference materials. SRM [967a](#) is a frozen human serum for evaluating the accuracy of procedures for the determination of creatinine in human serum. SRM [968e](#) is a lyophilized human serum for validating methods used to determine fat-soluble vitamins, carotenoids, and cholesterol in human serum and plasma. SRM [970](#) is a frozen human serum for validating methods for determining ascorbic acid in human serum and similar matrices. SRM [971](#) is a hormone in frozen human serum for evaluating the accuracy of procedures for the determination of the steroid hormones cortisol and progesterone in human serum. SRM [972](#) is a frozen human serum for evaluating the accuracy of clinical procedures for the determination of total cholesterol, HDL-cholesterol, LDL-cholesterol and triglycerides (triglycerides and total glyceride species). SRM [1950](#) is a frozen human serum for evaluating vitamin D metabolites. SRM [1951b](#) is a frozen human serum for evaluating the accuracy of clinical procedures for the determination of total cholesterol. SRM [1952a](#) is a freeze-dried human serum for evaluating PCB congeners, chlorinated pesticides, and PBDE congeners with non-certified values for PFCs and dioxins/furans. SRM [2921](#) is a human cardiac troponin complex. SRM 2921 is primarily intended for use in calibrating clinical procedures and devices for the determination of cardiac troponin I (cTnI) in human serum. It can also be used for value-assignment of calibrators and control materials.

For further information see: [SP 760.83](#)

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

SRM Description	909c	927d	956c	965b	967a	968e	970	971	972	1950	1951b	1952a
Unit Size	Human Serum (3 bottles x 2 mL each)	Bovine Serum Albumin (7% solution) (set (10))	Electrolytes in Frozen Human Serum (6 x 2.0 mL)	Glucose in Frozen Human Serum (set (8) (2 each conc))	Creatinine in Frozen Human Serum (set(4) (2 each conc))	Fat-Sol Vitamins, Carotenoids, and Cholesterol in Human Serum (set (3) (1 each conc.))	Ascorbic Acid in Frozen Human Serum (set (4) (2 each conc))	Hormones in Frozen Human Serum (2 x 5 mL)	Vitamin D in Human Serum (set (4) (1 each level))	Metabolites in Human Plasma (5 vials @ 1.0 mL each)	Lipids in Frozen Human Serum (set (4) (2 each conc))	Cholesterol in Human Serum (set (6) (2 each conc))
α-Carotene-Total						X				X		
α-Tocopherol						3 levels				X		
γ-Tocopherol						3 levels				X		
β-cryptoxanthin-total			X			3 levels						
Calcium Ionized			X									
5-Methyltetra-hydrofolic Acid										X		
Albumin		X										
Amino Acids										X		
Ascorbic Acid						2 levels						
Bilirubin										X		
Cardiac troponin C (cTnC)												
Cardiac troponin I (cTnI)												
Cardiac troponin T (cTnT)												
Chlorine Cl			X									
Cholesterol	X					3 levels				X	2 levels	3 levels
Cis-B-Carotene										X		
Coenzyme Q10						(3 levels)						
Cortisol								2 levels		X		
Creatinine	X				2 levels					X		
Dioxins/Furans												
Drugs of Abuse												
Elements (Selected)	X	3 levels								X		
Fatty Acids										X		
Folic Acid										X		
Glucose	X			4 levels						X		
Glycerides-Total	X									X	2 levels	
HDL-Cholesterol											2 levels	
Homocysteine										X		
LDL-Cholesterol											2 levels	
Lutein-Total						3 levels				X		
Lycopene-Total						3 levels				X		
Perfluorinated compounds (PFCs)										X		
Pesticides												
Polybrominated diphenyl ethers (PBDEs)												
Polychlorinated biphenyl (PCB) Congeners												
Progesterone								2 levels		X		
Prox. Protein	X	X								X		
Pyridoxal 5'-phosphate										X		
Retinol-Total						3 levels				X		
Retinyl palmitate										X		
Retinyl Stearate										X		
Selenium Species										X		
Testosterone								2 levels		X		
Total B-Cryptoxanthin										X		
trans-β-Carotene						3 levels				X		
trans-Lycopene						3 levels				X		
Triglycerides											2 levels	
Urea	X									X		
Uric Acid	X									X		
Vitamin D									4 levels			
Vitamin, 25-hydroxyvitamin D2										X		
Vitamin, 25-hydroxyvitamin D3						3 levels				X		
Vitamin, 3-epi-25-hydroxyvitamin D3												
Zeaxanthin-total						3 levels				X		

Certified values are normal font.
 Reference values are italicized.
 Values in parentheses are for information only.

These SRMs serve a variety of clinical measurement needs. SRM 900a is a frozen human serum for use in determining specific constituents. SRM 922a is a bovine serum albumin in a sterile 7% solution for use in the calibration and standardization of procedures to analyze total serum protein. SRM 456a is a frozen human serum for use in the calibration and standardization of procedures for the determination of specific electrolytes in either diluted or undiluted human serum. SRM 9065a is a frozen human serum for evaluating the accuracy of procedures used to determine glucose in human serum and to validate secondary reference materials. SRM 907a is a frozen human serum for evaluating the accuracy of procedures for the determination of creatinine in human serum. SRM 9684a is a lyophilized human serum for validating methods used to determine fat-soluble vitamins, carotenoids, and cholesterol in human serum and plasma. SRM 970a is a frozen human serum for validating methods for determining ascorbic acid in human serum and similar matrices. SRM 971a is a hormone in frozen human serum for evaluating the accuracy of procedures for the determination of the steroid hormones cortisol and progesterone in human serum. SRM 1971a is a frozen human serum for evaluating the accuracy of clinical procedures for the determination of total cholesterol, HDL-cholesterol, LDL-cholesterol and triglycerides (triglycerides and total glyceride species). SRM 277 is a frozen human serum for evaluating vitamin D metabolites. SRM 1971a is a frozen human serum for evaluating the accuracy of clinical procedures for the determination of total cholesterol. SRM 1558a is a freeze-dried human serum for evaluating PCB congeners, chlorinated pesticides and total cholesterol in human serum and similar matrices. SRMs 1485a and 1486a are freeze-dried human serums for evaluating PCB congeners, chlorinated pesticides, and PBDE congeners with non-certified values for PFCs and dioxins/furans. SRM 2921 is a human cardiac troponin complex. SRM 2921 is primarily intended for use in calibrating clinical procedures and devices for the determination of cardiac troponin I (cTnI) in human serum. It can also be used for value-assignment of calibrators and control materials.

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

[illegible]

Certified values are normal font.
Reference values are italicized.
Values in parentheses are for information only.