

### 106.3 - Mercury in Activated Carbon

These SRMs are intended for use in the evaluation of chemical methods of analysis for mercury.

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

SRM	2445	2448	2451
Description	Mercury in Iodized Activated Carbon	Mercury in Brominated Activated Carbon	Fine Carbon (Activated) - From Cyanide Ore Leaching
Unit of Issue	(25 g)	(25 g)	(100 g)
Mercury (mg/kg)	0.723	0.838	688

- Certified values are normal font
- Reference values are italicized
- Values in parentheses are for information only