

### 203.1 - Combustion Calorimetry (powder form)

These SRMs are for use as standards for calibration of combustion bomb calorimeters used in checking the performance of apparatus and analytical procedures.

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

| SRM  | Description  | Unit of Issue | Heat of Combustion (in MJ/kg)* |
|------|--|---------------|--------------------------------|
| 39j  | Benzoic Acid (Calorimetric Standard)                           | 30 g          | 26.434                         |
| 1656 | Thianthrene Combustion Calorimetric Standard                   | 30 g          | 33.480                         |
| 1657 | Synthetic Refuse-Derived Fuel Combustion Calorimetric Standard | 100 g         | 13.87**                        |
| 2151 | Nicotinic acid   | 25 g          | 22.184                         |
| 2152 | Urea Combustion Calorimetric Standard                          | 25 g          | 10.536                         |

\* The calorific values (MJ/kg) may decrease upon the aging or normal oxidation of the coals. NIST will continue to monitor these calorific values and report any substantive change to the purchaser.

\*\* Gross calorific value or HHV (Higher Heating Value).

- Certified values are normal font  
- Reference values are italicized  
- Values in parentheses are for information only