

3 - Thermal Analysis Purity Set

SRM 1514 is for evaluating methods of determining purity by differential scanning calorimetry. It consists of pure phenacetin and phenacetin doped with p-aminobenzoic acid.

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

3(3)- Thermal Analysis Purity Set

SRM	Description	Unit of Issue	Measurand
1514	Thermal Analysis Purity Set	set (4)	4 levels of p-ABA (in mol %)

- Certified values are normal font
- Reference values are italicized
- Values in parentheses are for information only