

106.3 - Mercury in Activated Carbon

This SRM is for use in the evaluation of chemical methods of analysis and in the calibration of instrumental methods of analysis.

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

SRM	Description	Unit of Issue	Mercury (mg/kg)
2445	Mercury in Iodinated Activated Carbon	25 g	0.743
2448	Mercury in Brominated Activated Carbon	25 g	0.838 mg/kg
2451	Fine Carbon (Activated) - From Cyanide Ore Leaching	100 g	688

- Certified values are normal font
- Reference values are italicized
- Values in parentheses are for information only