

106.6 - Used Auto Catalysts

These SRMs are intended for use in the evaluation of methods for the analysis of the platinum group metals and lead in auto catalysts. They were produced in cooperation with the International Precious Metals Institute and are issued as fine (74 μm , 200 mesh) powders.

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

Elemental Composition (in mg/kg)

SRM	Description	Unit of Issue	Lead (Pb)	Palladium (Pd)	Platinum (Pt)	Rhodium (Rh)
2556	Used Auto Catalyst (Pellets)	70 g	6228	326.0	697.4	51.2
2557	Used Auto Catalyst (Monolith)	70 g	13931	233.2	1131	135.1

- Certified values are normal font
- Reference values are italicized
- Values in parentheses are for information only