

201.2 - Biological Buffer Systems (powder form)

These SRMs are used to calibrate clinical instruments (e.g., blood pH measurements), in the physiologically important range of pH 7 to 8. They are based on a biological buffer system for clinical pH measurements and are each certified for use as an admixture only. The pH(S) values for the buffer solutions are certified at 0.05 molal and 0.08 molal with respect to the free acid and the sodium salt admixture as a function of temperature. The certified temperature range is from 0° C to 50° C.

PLEASE NOTE: The tables are presented to facilitate comparisons among a family of materials to help customers select the best SRM for their needs. For specific values and uncertainties, the certificate is the only official source.

SRM	Description	Unit Size	pH(S) Values (at 37 C)	
			0.05 molal	0.08 molal
2181	HEPES Free Acid	60 g	7.364	7.373
2182	NaHEPESate	60 g		
2183	MOPSO Free Acid	50 g	6.699	6.694
2184	NaMOPSOate	50 g		

Certified values are normal font.

Reference values are italicized.

Values in parentheses are for information only.